



True3 Digital Benefits

True Offloading through Sequential Whole-Body Pressure Relief

“Active” vs “Reactive” Technology - Alternating pressure relief regardless of patient’s movements or position.
Whole Body Sequential Pressure Reduction to ZERO every 7.5 minutes

TRI-PHASIC Alternation to promote Reactive Hyperemia

Sequential, “rolling” of 3 different pressures promotes capillary micro-bursts to help maximize blood flow.

Self-Contained Design – No Powerful Fan Blowing Air Constantly

Extremely Quiet for patient comfort and to support sleep - No added Heat to the room for patient comfort

Infection Control

No Bacteria/Virus spread from air sucked into pump or airflow out of cells, across patient, and into room

Fail Safe “3rd Row” of air cells Supports Patient even through Power Failures

Bottom rows of cells dedicated to patient protection via a dedicated fill tube with a one-way air flow valve.
Air cannot escape even if power is lost or pump is disconnected or malfunctioning.

3-Cell Support for More Comfort, Less Pain

4 inflated cells surround each deflated cell so that the patient is constantly supported.
Patient lies “On” the bed instead of being immersed “In” the bed for a more natural body alignment.

Staff Convenience

No Special Linens or Expensive Breathable Incontinence Pads Needed.
Simple Design – Strap Down Mattress, Attach Pump, Plug In & Turn On.
Prominent CPR Release at handle attachment to Pump

Reliability

FULL 2 YEAR WARRANTY against manufacturer defects